

# Ductless, Clean Bench

## Horizontal Clean Benches

Remove particles from room air to supply Class 100 air to your bench working area

- Provides ISO Class 5 conditions and product protection
- Two interior depths offered to match your application needs
- Easy to clean, static-dissipative stainless steel work surface

**NEW**

Sterile horizontal laminar flow benches filter out particles as small as 0.3 µm with 99.99% efficiency, providing ISO Class 5 air (fewer than 3520 particles, 0.5 µm or larger per cubic meter of air). Prefilters trap large particles, extending HEPA filter life. Differential pressure gauge displays pressure across the HEPA filter.

Benches are constructed of powder-coated steel exterior, safety glass side panels with utility ports, and 304 stainless steel work surface. Models with built-in 254 nm UV light feature retractable UV screen. Light timer provides timed UV light exposure (5, 10, 15, 30, 60, 120, 240 minutes and continuous settings). Benches include a built-in motor blower with vibration isolation pad. Other 230-volt plug types are available; contact an Application Specialist for details.

**Note:** Benches provide sample protection only. No personnel protection is provided. Do not use with biohazards, toxins or radionuclides.

**Compliance:** ETL (115 VAC models), CE (230 VAC models)

**What's included:** 10-ft power cord with plug.



Bench 09114-12 shown with base stand 09116-76

ISO 9001:2008  
CERTIFIED SUPPLIER



Freight  
Transport  
Items over 100 lb

Nominal width (ft)	Dimensions (W x H x D) <sup>†</sup>		Power VAC, Hz	Plug	Shpg wt lb (kg)	Clean benches		Clean benches with UV light	
	in.	cm				Cat. no.	Price	Cat. no.	Price
3	38 x 48.4 x 34	96.5 x 123 x 86.4	115, 60	115 VAC, 15 A	364	T-09116-00		T-09116-32	
			230, 50/60	British (UK)	(165)	T-09116-02		T-09116-34	
3	38 x 48.4 x 39	96.5 x 123 x 99	115, 60	115 VAC, 15 A	372	T-09116-04		T-09116-36	
			230, 50/60	British (UK)	(169)	T-09116-06		T-09116-38	
4	50 x 48.4 x 34	127 x 123 x 86.4	115, 60	115 VAC, 15 A	455	T-09116-08		T-09116-40	
			230, 50/60	British (UK)	(206)	T-09116-10		T-09116-42	
4	50 x 48.4 x 39	96.5 x 123 x 99	115, 60	115 VAC, 15 A	465	T-09116-12		T-09116-44	
			230, 50/60	British (UK)	(211)	T-09116-14		T-09116-46	
6	74 x 48.4 x 34	188 x 123 x 86.4	115, 60	115 VAC, 15 A	645	T-09116-16		T-09116-48	
			230, 50/60	British (UK)	(293)	T-09116-18		T-09116-50	
6	74 x 48.4 x 39	188 x 123 x 99	115, 60	115 VAC, 15 A	660	T-09116-20		T-09116-52	
			230, 50/60	British (UK)	(299)	T-09116-22		T-09116-54	
8	99 x 48.4 x 34	251.5 x 123 x 86.4	115, 60	115 VAC, 15 A	835	T-09116-24		T-09116-56	
			230, 50/60	British (UK)	(379)	T-09116-26		T-09116-58	
8	99 x 48.4 x 39	251.5 x 123 x 99	115, 60	115 VAC, 15 A	855	T-09116-28		T-09116-60	
			230, 50/60	British (UK)	(388)	T-09116-30		T-09116-62	

<sup>†</sup>For benches with UV light, add 1" (2.5 cm) to the width.



09116-80

## Replacement Filters

Bench width	Prefilters		HEPA filters	
	Cat. no.	Price	Cat. no.	Price
3 ft	T-09114-90 <sup>‡</sup>		T-09116-82	
4 ft	T-09114-91 <sup>‡</sup>		T-09116-84	
6 ft	T-09114-91 <sup>‡</sup>		T-09116-86	
8 ft	T-09114-91 <sup>‡</sup>		T-09116-84	

<sup>‡</sup>The 38 to 39" W and 50 to 51" W clean benches require two prefilters, 74 to 75" wide clean benches require three prefilters, and 99 to 100" wide clean benches require four prefilters.

<sup>††</sup>The 99-100" wide clean benches require two HEPA filters.

## Labconco Telescoping Base Stands

Stands feature telescoping legs so you can select stand height during installation. Height is adjustable to eight positions in 1" increments from 27.5" to 34.5". Frame is durable 1.75" tubular, powder-coated steel with rear crossbar supports. ADA compliant. Fixed leveling feet. 1000 lb (454 kg) load capacity.

Catalog number	For nominal width	Base stand dimensions (W x D)		Shpg wt lb (kg)	Price
		in.	cm		
T-09116-74	3 ft	35 x 29	89 x 73.7	100 (45)	
T-09116-76	4 ft	47 x 29	119.4 x 73.7	120 (55)	
T-09116-78	6 ft	71 x 29	180.3 x 73.7	135 (61)	
T-09116-80	8 ft	96 x 29	243.8 x 73.7	165 (75)	